



Installation Guide Executive Document Note: Study iPhone password: 123456, Apple Watch: 1234

Science CHING FURTHER

Pairing the Apple Watch

- Press and hold the side button until you see the Apple logo -
- Bring the iPhone near the Apple Watch, wait for the Apple Watch pairing screen to appear on your iPhone then tap "continue"
- Re-open the Apple Watch App on your iPhone, then tap "Pair New Watch". _



- When prompted, position your iPhone camera over the Apple Watch appears so that it appears _ in the viewfinder of the Apple Watch App. This pairs the two devices.
- Tap "Set Up Apple Watch". Follow the instructions on your iPhone and Apple Watch to finish setup.



Choosing a Watch Face

- Open "Watch" App on the iPhone
- On the "Face Gallery" tab, scroll down to "Infograph Modular" watch faces _







Installation Guide Executive Document



- Choose from the following complications, then click "Add"
- There should be 2 complications (all the rest should be "off")
 - Top left: Battery
 - Bottom middle: Back2Play

*This will be "off" if the B2P App is not yet installed

- On the "My Watch" tab, ensure that the Watch face you just created is the only option under "My Faces"
- If there are multiple watch face options, click "edit" then use the minus sign button to remove all watch faces except those selected (Battery and B2Play App)

Face Gallery	Û
······································	Infograph Modular
() ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	
Complications	
Top Left	Battery >
Date	Off >
Middle	Off >
Bottom Left	Off >
Bottom Middle	Back2Play >
Bottom Right	Off >
O Free	2

Personal iPhone setup

- Ensure you have an iPhone 8 or later that has been updated to iOS 13 or later
- Ensure you have an Apple ID that has been logged in







Installation Guide Executive Document

Back2Play App download on the Apple Watch

- Ensure the Apple Watch has been paired to the iPhone and the iPhone is running on iOS 13 or later BEFORE you download the Back2Play App
- Ensure "My Watch" App is on the iPhone
- On "My Watch App" scroll until you see the "Back2Play App", then click "Install"
- Add the Back2Play App watch complication
- Open "My Watch" app from the iPhone
- On the "My Watch" tab, click the Watch face you created
- For the "Bottom Middle" select Back2Play
- Open the App on the phone and enter the Back2Play App access code
- Make sure you allow the Back2Play App access to Health data, notifications and motion/ accelerometer
- Check that the accelerometer has permission to track movement. On the phone go to Settings>Privacy>Motion and Fitness and make sure Fitness Tracking, Back2Play and Health are turned on
- Open the App on the Watch, this is required so accelerometer data is tracked.